ENDOCLOT™

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
# EndoClot™ Polysaccharide Hemostatic System

**Caution:** EndoClot™ Polysaccharide Hemostatic System should only be used by physicians or other licensed practitioners.

## DESCRIPTION

EndoClot™ Polysaccharide Hemostatic System (EndoClot™ PHS) is a medical device composed of absorbable modified polymer (AMP®) particles and an applicator with a tube connecting to a clean external gas source. AMP® particles are biocompatible, non-pyrogenic and derived from plant starch. The device contains no human or animal components. EndoClot™ PHS is intended as an adjunct hemostat to control bleeding from capillary, venous, or arteriolar vessels in endoscopic procedures when conventional procedures are ineffective or impractical.

**Schematic diagrams of AMP® particles and Applicator**

## ACTION

AMP® particles have a molecular structure that rapidly absorbs water from the blood. This dehydration process causes a high concentration of platelets, red blood cells, and coagulation proteins (thrombin, fibrinogen, etc.), which accelerates the normal, physiologic clotting cascade. When in contact with blood, AMP® particles support the formation of a gelled, adhesive matrix which provides a mechanical barrier to further bleeding. AMP® particles are degraded within an hour to several days depending on the amount of material applied and the site where it is used. AMP® particles are degraded by amylase and glucoamylase.

## INDICATIONS

EndoClot™ PHS is indicated for use in endoscopic procedures as an adjunct hemostat when the control of bleeding from capillary, venous or arteriolar vessels by conventional procedures is ineffective or impractical.

## INSTRUCTIONS FOR USE

The following instructions do not eliminate the necessity of formal training for the use of EndoClot™ PHS. The techniques and procedures described here do not represent all clinical acceptable protocols, nor are they intended as a substitute for a physician's experience and judgment in treating specific surgical conditions.


Fig. 1


Fig. 2


Fig. 3


Fig. 4

1. Screw the AMP® particle dispenser firmly onto the end of the gas/powder mixing chamber (Fig.1).
2. Connect the gas filter to an air source with constant airflow ( $18 \pm 3$  kPa) to prevent fluid from reverse flowing into the catheter (Fig.2).

3. Insert the catheter through a working channel on the endoscope (Fig.3) and position the endoscope to the bleeding site. Avoid direct contact of the catheter tip with any fluids.
4. Select a gas source with constant airflow ( $12 \pm 3$  kPa). Gently tilt the dispenser to deliver the powder to the site of bleeding via the gas/powder mixing chamber and catheter (Fig.4).
5. Remove the catheter from the endoscope immediately after hemostasis is achieved.

### **Application Techniques**

1. Identify the bleeding lesion(s). Removing excess blood from the site of bleeding is essential to achieve optimal hemostatic efficacy.
2. Insert the catheter into the endoscope and position its tip to the site of bleeding. Spray a liberal amount of AMP® particles on the wound.
3. Direct contact of the catheter tip with any fluid must be avoided to prevent blockage caused by reverse flow. If the catheter is occluded, replace it with a new one then repeat the process.
4. If bleeding continues, irrigate excess AMP® particles with a saline solution and aspiration, then re-apply.
5. Once hemostasis is achieved, remove excess AMP® particles with irrigation and aspiration.
6. Following the procedure, extract gas from gastrointestinal tract through the endoscope and ensure the endoscope is completely cleaned by irrigation to avoid endoscope working channel occlusion.
7. Hold the dispenser and gas/powder mixing chamber in an upright position to avoid spillage when inserting the catheter into the endoscope.
8. Avoid folding the catheter when delivering the powder.

### **CONTRAINDICATIONS**

EndoClot™ PHS is contraindicated in patients who are sensitive to starch or starch-derived materials.

Do not inject AMP® particles into blood vessels: extensive intravascular coagulation may occur.

Do not inject AMP® particles into the bile duct or pancreatic duct.

Do not inject AMP® particles into the bladder or ureteral lumen.

Do not use AMP® particles for controlling post-partum bleeding or menorrhagia.

## WARNINGS

EndoClot™ PHS is not intended as a substitute for good surgical practice, and in particular, the proper use of conventional procedures for hemostasis.

Combined use of EndoClot™ PHS with other topical hemostatic agents has not been studied in controlled clinical trials.

EndoClot™ PHS is for single use only.

Apply quickly; deflation is required when EndoClot™ PHS is used in the gastrointestinal tract for over a minute.

## PRECAUTIONS

EndoClot™ PHS is not recommended as a primary treatment for coagulation disorders.

EndoClot™ PHS is intended to be used in a dry state. Contact with fluids prior to application will result in the loss of hemostatic efficacy.

## ADVERSE REACTIONS

None reported to date.

## DOSAGE AND ADMINISTRATION

Aseptic technique should always be used. A liberal amount of AMP® particles should be applied to the bleeding site until hemostasis is achieved. After hemostasis is achieved, AMP® particles should be removed by irrigation and/or aspiration.

## HOW SUPPLIED

EndoClot™ PHS			
Ref No.	AMP® Particles	Catheter Length	Packaging
EPK1801	1g	1800mm	1 set/box
EPK1802	2g	1800mm	1 set/box
EPK1803	3g	1800mm	1 set/box
EPK2301	1g	2300mm	1 set/box
EPK2302	2g	2300mm	1 set/box
EPK2303	3g	2300mm	1 set/box
EndoClot™ Catheter			
EPAA180	-	1800mm	5 pcs/box
EPAA230	-	2300mm	5 pcs/box

## STERILIZING METHOD & EXPIRATION DATE

EndoClot™ PHS is sterilized and should not be re-sterilized. The AMP® particles are sterilized by Irradiation. EndoClot™ applicator is sterilized by EO. Discard the unused, open packages properly.

If stored under the conditions specified in this manual (see STORAGE AND HANDLING), the unopened and undamaged product remains sterile for three years from the date of sterilization.

## STORAGE AND HANDLING

EndoClot™ PHS should be stored under room temperature. EndoClot™ PHS should not be exposed to temperature ranges outside -40°C to 60°C during transportation. EndoClot™ PHS should be used immediately after the package has been opened.

## DISPOSAL

This product should be disposed of in compliance with pertinent government regulations regarding medical devices.

## LIMITED WARRANTY

EndoClot Plus, Inc. warrants that this product is free from defects in workmanship and materials. Liability under this warranty is limited to refund or replacement of any product, which has been found by EndoClot Plus, Inc. to be defective in workmanship and materials. EndoClot Plus, Inc. shall not be liable for damages arising from the use, misuse, or abuse of this product or its content in ways that are inconsistent with the specific indications described in these Instructions for Use. Damage to the product through misuse, alteration, improper storage, or improper handling shall void this limited warranty.

No employee, agent, or distributor of EndoClot Plus, Inc. has authority to alter this limited warranty in any respect. Any purported alteration or amendment shall not be enforceable against EndoClot Plus, Inc., and should be reported to EndoClot Plus, Inc. and/or appropriate authorities.

THIS WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER OBLIGATION OF ENDOCLOT PLUS, INC.

## TRADEMARKS

EndoClot™ is a trademark of EndoClot Plus, Inc.

AMP® is a trademark of Starch Medical Inc.

= Do not reuse

= Do not resterilize

= Do not use if package is damaged

= Use-by date

= Catalogue number

= Sterilized using irradiation

= Sterilized using ethylene oxide

= Batch code

= Date of manufacture

= Caution

= CE-mark and identification number of Notified Body.

Certified according to MDD (93/42/EEC)

= Authorized representative in the EC

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